

a1 medical



CTS

CORNEA TRANSPLANT SYSTEM

KRUMEICH TREPHINE



C-1100

CTS Vorderkammerbank-System, komplett
CTS Anterior Chamber System, complete



C-1200

CTS Trepanations-System
CTS Trepanation-System




C-1300

CTS Vakuum-Saugring
CTS Vacuum-Suction-Ring

**C-1100**


CTS Vorderkammerbank-System, komplett, bestehend aus
CTS Anterior Chamber System, complete, consisting of

**C-1107**

CTS Spannschraube für Vorderkammerbank
CTS Tensioning Screw for Anterior Chamber

**C-1105**

CTS Fixier-Ringplatte, Innen-Ø 11.0 mm
CTS Fixation Ring Plate, inner-Ø 11.0 mm

**C-1106**

CTS Fixier-Ringplatte, Innen-Ø 12.0 mm
CTS Fixation Ring Plate, inner-Ø 12.0 mm

**C-1102**

CTS Zentralpils mit O-Ring
CTS Central Plug with o-ring

**C-1101**

CTS Vorderkammerbank, allein
CTS Anterior Chamber, only



C-1200

CTS Trepanations-System, komplett, bestehend aus
CTS Trepanation-System, complete, consisting of



C-1202

CTS Trepanker
CTS Trephine Core



C-1201

CTS Trepanations-Hülse für Trepanker
CTS Trepanation-Guide for the Trephine



Trepan
Trephine

C-1270 Ø 7.0 mm

C-1275 Ø 7.5 mm

C-1280 Ø 8.0 mm



Glas-Obturatoren-Set
Set of Glass-Obturators

C-1271

C-1276

C-1281

mit folgenden Radien und Diopter
with the following radii and diopter

R: 7.0 mm	R: 7.1 mm	R: 7.2 mm	R: 7.3 mm
D: 47.56	D: 46.89	D: 46.24	D: 45.60
C-1271+7.0	C-1271+7.1	C-1271+7.2	C-1271+7.3
C-1276+7.0	C-1276+7.1	C-1276+7.2	C-1276+7.3
C-1281+7.0	C-1281+7.1	C-1281+7.2	C-1281+7.3

R: 7.4 mm	R: 7.5 mm	R: 7.6 mm	R: 7.7 mm
D: 44.99	D: 44.39	D: 43.80	D: 43.23
C-1271+7.4	C-1271+7.5	C-1271+7.6	C-1271+7.7
C-1276+7.4	C-1276+7.5	C-1276+7.6	C-1276+7.7
C-1281+7.4	C-1281+7.5	C-1281+7.6	C-1281+7.7

R: 7.8 mm	R: 7.9 mm	R: 8.0 mm
D: 42.68	D: 42.14	D: 41.61
C-1271+7.8	C-1271+7.9	C-1271+8.0
C-1276+7.8	C-1276+7.9	C-1276+8.0
C-1281+7.8	C-1281+7.9	C-1281+8.0

C-1300

CTS Vakuum-Saugring, komplett, bestehen aus
CTS Vacuum-Suction-Ring, complete, consisting of



C-1301

CTS Saugring
CTS Suction-Ring



C-1302

CTS Vakuum-Ring-Einsatz
CTS Vacuum-Ring-Insert

Zusätzliche Bestandteile des CTS – Cornea Transplant Systems / Krumeich Trepan

Additional Components of the CTS – Cornea Transplant System / Krumeich Trephine



C-1303

CTS Zentrier-Ring
CTS Centering Ring



C-1108

CTS DMEK-Teflon-Pilz
CTS DMEK-Teflon-Central Plug



C-1401

CTS Trepan, Ø 15.0 mm, zum Vortrepanieren der Spederhornhaut
CTS Trephine, Ø 15.0 mm, for the Pre-Trepanation of the Donor Cornea



C-1400

CTS Schneidblock zum Vortrepanieren der Spederhornhaut, Teflon
CTS Cutting-Block for the Pre-Trepanation of the Donor Cornea, Teflon



MIS-2000

CTS Sterilisationskassette, komplett, Abmessungen: 277x181x50 mm
CTS Sterilization Case, complete, Dimensions: 277x181x50 mm



C-1104

CTS T-Schraubendreher für Zentralpilz und DMEK-Teflon-Pilz
CTS T-screwdriver for Central Plug and DMEK-Teflon Plug



C-1500

CTS Tonometer, 22.0 mm HG
CTS Tonometer, 22.0 mm HG



F-4135

CTS Pinzette zum Einsetzen und Entfernen des Vakuum-Ring-Einsatzes
CTS Forceps for inserting and removing the Vacuum-Ring-Insert

